



Acoustic (Non-contact) Liquid Dispenser



- Affordable non-contact, acoustic liquid dispenser. Low energy acoustics used, safe for all biological fluids (including cells).
- Capable of plate reformatting, compound library screening, assay development, cherry picking and serial dilutions (1C50s) in and out of 384, 1536 and 3456 well plates.*
- Precise and controlled transfer of chemical and biological compounds. Suitable for dispensing wide range of fluids e.g. 100% DMSO, aqueous buffers, glycerol solutions. Standard unit able to dispense between 1.0 - 250nl.
- The Ats -100 creates non-contact quantitative microarrays, even on top of existing arrays. Suitable for use with proteins, tissues and whole blood. This option is able to dispense volumes between 250pl - 250nl.
- Can be easily integrated to produce a complete high throughput screening system.
- Low dead volumes. e.g.<1.0ul for 1536 wellplate.
- Protocols are remotely configured and downloaded to the Ats-100 for processing. Software is designed to handle multiple source plates, only one file required for cherry picking.
- Continuously drain wells from top to bottom, independent of neighbouring liquid levels. Real-time measurement of all liquid levels in all wells performed.

* Certified for use with COC plates from Aurora Discovery and Greiner Bio-One. Contact Cedrex for full list of compatible plates.

Cedrex A/S
Notgangen 3
DK-2690 Karlslunde
Phone +45 3257 5181
Fax +45 3257 5182
Info@cedrex.dk
www.cedrex.dk



Acoustic (Non-contact) Liquid Dispenser

SPECIFICATIONS

Fluid Viscosity	0.05-5 cP
Surface Tension	27-72 dynes/cm
Density	0.5-2.0 g/cm ³
Transfer Volume	1-250nl 250 pl - 250nl with microarray package
Transfer Variability	<10%CV
384 well plate replication 1536 well plate replication 3456 well plate replication	<2min* <6min* <15min* *for single unit volume transfer per well
Dimensions	W=700m x D=760mm x H=700mm
Weight	105kg
Electrical Requirements	240Vac, 50Hz
Clean Dry Air	60-80 psi, 2CFM
System Components	Internal coupling fluid circulation system Well plate vision system Windows-based PC Internal bar code reader Internal target plate flipper and deionizer Ethernet/connectivity data Real-time data collection and management system

Ats-100 manufactured by EDC Biosystems
www.edcbiosystems.com